

## AUTOMATIC BOAR SEMEN COLLECTION



Available with electrical lifting columns with control units

**FUTURO** is a new product for automatic boar semen collection, that allows to reduce the time needed for collection, while enhancing its efficiency thanks to an automated system suitable for all kinds of boars. It is easy to use and can be fit into existing equipment respecting the highest hygienic standards during the collection.

The height of the collection dummy can be adjusted by means of a convenient electric piston. The boar jumps on it in rest position (the lowest one), then the operator adjusts the height according to the dimensions of the boar itself by simply pressing a button.

The boar will find the ideal position and comfort to ensure abundant ejaculates.

Once collection is over, the collection dummy goes back to its rest position (the lowest one), thus facilitating the descent of the boar.

The electric adjustment of the collection dummy height makes the automatic collection system suitable for all kinds of boars, being them young or adult.

**FUTURO's** main component is the artificial texture cervix, that perfectly simulates the natural one.

The sheath material ensures a good and comfortable catch for the boar and facilitates the initial collection stages, during which the operator puts the artificial cervix in place.

The sheath thickness and softness prevent any injury to the boar penis.



1. Sliding System
2. Rubber Pad o Collection Sheath
3. Sheath Forceps
4. Themos Restraint Device
5. Stainless Stell Dummy
6. Height Adjustment
7. Support

The rubber pads gently fit into the boar penis, ensuring the right catch. The **FUTURO** collection dummy is supplied with a sliding system enabling the boar to move forwards and backwards during the collection.

After prepuce cleaning, at the

beginning of the erection, the operator manually discards the first part of the collected material; immediately after, he fixes the penis to the artificial cervix by means of a special fork, connecting at the same time the collecting vat, which is equipped with a filter-provided bag.

The collection can then proceed automatically until its conclusion, allowing the operator to deal with other activities such as the preparation of the next boar.

The **FUTURO** collection system does not require any maintenance.

In designing the automatic collection system, the utmost importance was given to user-friendliness and limited maintenance requirements, as well as to the efficiency of the system over time.

### Advantages

- Time saving up to 80% (1 supervisor every 3 collection dummies)
- Safety and comfort for the operator
- Very easy collection, even for inexperienced operators
- Abundant ejaculate
- Utmost comfort and safety also for the boars
- Lower labour costs per produced dose
- Suitable for all kinds of boars (both young and adult ones)
- Reduction of bacterial contamination



Available in 2 Versions: Plastic Coated and Stainless Steel

### FUTURO COLLECTION SYSTEM

Code	Description
P60300	FUTURO BASE stainless stell
P60400	FUTURO COMPLETE stainless stell
P60310	COLLECTION SHEATH
P60350	SHEATH FORCEPS
P60311	SLIDING SYSTEM
P60312	Electrical lifting columns with control units
P60313	Manually lifting columns