




Sheathfix Universal Syringe 0,25 + 0,5

Universally suitable for all kinds of straws and sheaths.
Comfortable working through ergonomic design
AI instrument with best available length: only 1 cm of the sheath length is required for fixation
Simple cleaning: the handle is made of only one piece.
Sterilization in the autoclave (up to 180°C)
No loose parts, the plunger is secured against sliding out

SHEATHFIX UNIVERSAL SYRINGE

Code	Description
BM000025	Sheathfix Universal Syringe 0,25 + 0,5



UNIVERSAL MCI SHEATHS 50/bag, 1000 box

Suitable for all types of semen straws
Absolutely tight fit with all types of straws
Gentle insemination due to a soft shape of the tip.

MCI SHEATHS

Code	Description
BM000025	Sheaths 50/bag, 1000 box




NEW Thermostatic unit for BOVINE AI Syringes

Avoid thermal shocks that could prevent an optimum AI results It ensures a constant temperature of 37°.
The system includes 2 main independent components:
- Unit control 12V with 1 PID temperature regulators with digital Display.
- Container whit rack for 10 or 20 AI Syringes

THERMOSTATIC UNIT AI SYRINGES

Code	Description
BMOP10	Thermostatic unit for BOVINE AI Syringes



Thawing Unit MCI 30/60

- Electronic temperature control
- Available whit 30 mm or 60 mm diameter chamber
- Central, built-in thawing chamber with lift
- Preset temperature setting: +37°C,
- Accuracy: ± 0.2°C
- Power supply: 230 V
- Power supply: 110 V, 12 V DC (car battery) optional

THAWING UNIT MCI 30

Code	Description
BM000030	Thawing Unit MCI 30
BM000060	Thawing Unit MCI 60



MediChimica STRAW










Very accurate production techniques and our strict quality control make sure that your demands for a high standard will be consequently realised.













Medi Chimica straws are available in transparent colours.

The benefits is more sperm cells per ejaculate available

MEDICHIMICA STRAW 0.25 cc

Code	Description
 BMS000W25	White
 BMS000T25	Transparent
 BMS000R25	Red
 BMS000P25	Pink
 BMS000G25	Green
 BMS000B25	Blue
 BMS000U25	Turquoise
 BMS000P25	Purple
 BMS000G25	Grey
 BMS000Y25	Yellow
 BMS000O25	Orange
 BMS000S25	Salmon


MEDICHIMICA STRAW 0.50 cc

Code	Description
 BMS000W50	White
 BMS000T50	Transparent
 BMS000R50	Red
 BMS000P50	Pink
 BMS000G50	Green
 BMS000B50	Blue
 BMS000U50	Turquoise
 BMS000P50	Purple
 BMS000P50	Grey
 BMS000Y50	Yellow
 BMS000O50	Orange
 BMS000S50	Salmon



Individually folded for pop-up dispensing

Material : PE



SUPERSENSITIVE

A "silky" copolymer glove offering excellent tactile sensitivity. Ideal when precise touch sensitivity is essential.

SENSITIVE

Specially developed for artificial insemination, this polyethylene protection glove offers very good tactile sensitivity to make easy work of gestation checks.

ECONOMICAL

The most economical of our polyethylene gloves for veterinary purposes

DISPOSABLE GLOVES

ref.	Code	Description
1	P80100	SUPERSENSITIVE gloves 100 pieces package
2	P80101	SENSITIVE gloves 100 pieces package
3	P80102	ECONOMICAL gloves 100 pieces package



Medi Chimica β -Ketone

The Medi Chimica β -Ketone sensors can read blood β -Ketone levels between 0.0 and 8.0 mmol/L.

The blood β -Ketone test measures Beta-Hydroxybutyrate (β -OHB), the most important of the three ketone bodies in the blood .

Normally, levels of β -OHB are expected to be less than 0.6 mmol/L. Level up to 1,0 mmol/L detect a subclinical ketosi in bovine.

Advantages:

- HI analytical accuracy
- Rapid and easy test on field
- Only 10 second for analysis
- HI correlation and accuracy whit laboratory test

Medi Chimica β -Ketone

Code	Description
BMK0000	Medi Chimica β -Ketone