

MEDISCANMini

Digital Veterinary Scanner



Digital MEDISCAN-Mini, is the new contribution by MEDICHIMICA INTERNATIONAL for veterinary. It has the transducer, display and power supply integrated in one compact sturdy housing. The total unit weights 500g (including the battery) only. After use, it can be simply be cleaned under running tap water.

System Feature

- Digital technology
- 3.5 inch TFT color LCD
- Completely washable
- Simple to operate
- All in one, no cable
- 3hrs continuous operation
- 2 years warranty

Physical Characteristic

- Size: 171x93x62 mm
- Weight: 500g , including the battery

Standard Accessories

- AC-adapter (110/220V)
- two rechargeable battery
- One AA battery holder
- One battery charging cable
- One carrying bag
- One bottle of ultra gel
- One hand strap

Specifications

- Monitor: digital 3.5 inch TFT color LCD
Probe: Multi-frequency mechanical sector probe
Scan Frequency: 2.5MHz, 3.5MHz, 5.0MHz
Display Depth : 160mm
Scan Angle: 78 degree
Gray Shade: 256 levels
Measurement: Distance, Area, Volume, Heart Rate
OB Measurement: Swine-HL, Equine, Bovine-BL,SL,HL, Sheep-USD
Cat-HD,BD, Dog-BD, GSD, CRL,HD, BD
Permanent Storage: 127 frames
- AC-adapter: 110V-240V, 50/60Hz
Battery: (1 800mAh) Internal rechargeable battery
More than 3hrs for continuous operation
less than 3hrs for recharging

