

MEDISCAN HI

Palmsmart best image quality

weigh only 1.8kg



Powerfull function and original design

256 gray scales, 10,4 inches TFT LCD, internal battery.

High integrated circuits are applied for this type system, so it is small and easily carried.

Full alphanumeric, backlit keyboard.

From horses, cattle, sheep, goats and swine, high image quality and a wide range of available probes make the MEDISCAN HI an asset to any modern veterinary practice.

MEDISCAN HI

scan your future

MEDI
chimica
INTERNATIONAL



Palmsmart MEDISCAN HI is the latest model by MEDI CHIMICA INTERNATIONAL, weighting only 1.8kg, is distinguishable by its bright and clear image quality.

The interchangeable probes offer a wide range of diagnostic applications.

This scanner is a cost effective solution for family doctors, veterinary scanning.

The complete software and the high product versatility are available in it.



Model : MEDISCAN HI

- Scanning Method:Conves,linear triple-frequency with wide band.
- Display Mode:B,B+B,B+M,M
- Gray Scale:256 levels
- Monitor:10.4" TFT color LCD
- Image Procession:DFS,DRF,RDA,VGA/DSC and Image Post Procession
- 206 frames for cine-memory
- Image magnification in four levels: $\times 1.0, \times 1.2, \times 1.5, \times 2.0$ times
- Multi-step Zoom (1:2) is available in both real-time and frozen modes
- 4 focusing combination selectable
- Depth Shift is available in read-time in B and B+B mode
- Image Conversion is U/D, L/R and B/W
- 8 difference coded-color available on a external color monitor
- 15 difference kinds of body marks with defference calipers.
- Measurement is Distance, Circumference/Area, Time/Heart Rate EF rate for heart,and
- GS,CRL,BPD,HC,FL,AC for OB/GY measurement, and Gestational Age,Fetal weight,EDD calculation.
- Characteristic showing in real-time clock, ID, Sex, Age, Focusing info and Hospital Name

Medi Chimica International has an international distribution network. Contact us to have the address and phone number of your local distributor or if you want to become a distributor.

Medi Chimica International ensures you an efficient after sales service.

Medi Chimica International also provides training for operators who can learn how to interpret ultrasound scanning image and how to handle the probe on different species.