

# MEDISCANPRO with handle and internal battery

Portable, battery powered hand held ultrasound system with detachable probe heads for early Preg checking

*Easy to learn, easy to use and lightweight*



Real dimension

New probe with freeze button

## Lightweight

600 gr mainframe with small internal battery (1.30 lbs) allows for easy moving around the animals. Because the MEDISCAN PRO is extremely compact you can screen between 150 and 200 sows an hour

## Easy to operate

Logical and easy to reach controls make for easy operation and simple to learn to use

## Display quality

Built-in high resolution TFT LCD display offering high quality images and excellent contrast quality.

# MEDISCAN PRO

scan your future

**MEDI**  
**chimica**  
INTERNATIONAL



**ONLY 600 gr (1.3 lbs)**

- Early Preg Checking starting from 19-20 days (lower your non-productive days)
- Checking Fetal viability
- Identifying retained piglets
- Easy loin depth measurement
- Accurate backfat measurement

#### NEW PROBE - FEATURES

- Fully digital interface
- High frame rate
- High resolution
- Multi-frequency high resolution transducer
- Low noise
- Lower power consumption
- Excellent penetration
- Highly flexible cable



Also available MEDISCAN  
Version with external battery



It is the smallest B mode ultrasound diagnose system thanks to its high integrated circuits.

Adopting TFT 5 inch LCD and it is very low energy consumption.

With an exceptional quality of image fully adapted to pregnancy diagnosis, MEDISCAN PRO offers you unbeatable comfort in work for great efficiency.

#### Model : MEDISCAN PRO

- Standard probe: 3.5MHz mechanical sector probe;
- Scanning depth: 192 mm;
- Monitor: TFT 5.0 inch LCD;
- Battery: 2200mAh ,continuous working 2+hours;
- Display modes: B, B+B, B+M, M;
- Gray scales: 256;
- Image magnification:  $\times 1.0$ ,  $\times 1.2$ ,  $\times 1.5$ ,  $\times 2.0$ ;
- Measurement: Distance, Circumference, Area, Heart rate & Pregnant week;
- Characters: ID ,age ,sex;
- Power consumption: 13VA;
- Power supply: 10-12V battery / AC-adapter;
- Power adapter: 110/220 VAC system power adapter and battery charger;
- Weight: 600 grams (1.3 lbs);
- Dimensions: 23,6 x 4,1 x 12,9 cm ( $9.2 \times 1.65 \times 4.8$  inches).

Medi Chimica International has an international distribution network. Contact us to have the address and phone number of your local distributor or if you want to become a distributor.

Medi Chimica International ensures you an efficient after sales service.

Medi Chimica International also provides training for operators who can learn how to interpret ultrasound scanning image and how to handle the probe on different species.