


MEDI Needle

**STRENGTH, SHARPNESS,
NEW PACKAGING
EASY IDENTIFICATION**



One-at-a-time pop-pup dispensing



Medi Needle (Steel)



CE

Medi One (Plastic)

MEDI NEEDLES CHARACTERISTICS (ALUMINIUM HUB)

- All aluminium manufactured
- 100 pcs per box
- Single packaging
- Disposable needles
- Size stamp in each needle

STRENGTH OF THE MEDI NEEDLE LINE:

Strength:

Needles are provided with thick-walled cannulas, which are crimped and pressed into sturdy plated brass hubs. These stronger sidewalls increase load strength. The bending and breaking risks are six times lower than with conventional needles, and needle tips remain sharp for a longer time. These superior features make these needles last three times longer than conventional injection needles.

Sharpness:

Needles are ultra-sharp, with three-point precision ground tips for sharper cuts and easier penetration. Sharpness increases the number of possible injections per needle and reduces the risk of injection site abscesses.

Packaging:

Innovative transparent container suitable for carriage and deposit. Needles are packed in colour-coded cartridges for safe and easy transport. Cartridge labels are tamper evident and identify needle length and gauge.

- Detectable thick wall stainless steel cannula available

MEDI NEEDLES CHARACTERISTICS (STAINLESS STEEL)

- One-at-a-time pop-pup dispensing
- Needle colour code for quick and easy identification of the diameter
- Innovating transparent carrying and deposit container
- Each plastic box contain 12 needles and stainless steel bar to clean
- 120 Needles per carton box
- Sterilizing needles
- Dimension on the label are in Gauge and mm
- Manufactured from a springy drawn nickel chromium stainless steel
- Thick walled needle tube with a very high stability and elasticity
- Smooth surface tube
- Needle mount with a generous styled supporting surface
- Tube diameter stamped on needle mount
- Distinct bison bevel



CE

Medi One (Aluminium)

MEDI NEEDLES CHARACTERISTICS (PLASTIC HUB)

- Plastic hub Dimension on the label are in Gauge and mm and metal needle
- 100 Pcs per box
- Plastic hub with different colour
- Colour code for quick and easy identification of the diameter
- Dimension on the label are in Gauge and mm
- Disposable needles

MEDI Needle

STRENGTH, SHARPNESS,
NEW PACKAGING
EASY IDENTIFICATION

MEDI Needle



MEDI NEEDLE (STEEL)

Code	Lock	Sizes
MPLL1010	Luer Lock	1,0 x 10 19 G x 3/8" *
MPLL1015	Luer Lock	1,0 x 15 19G x 5/8" *
MPLL1020	Luer Lock	1,0 x 20 19G x 3/4" *
MPLL1210	Luer Lock	1,2 x 10 18G x 3/8" *
MPLL1215	Luer Lock	1,2 x 15 18G x 5/8" *
MPLL1220	Luer Lock	1,2 x 20 18G x 3/4" *
MPLL1225	Luer Lock	1,2 x 25 18G x 1" *
MPLL1415	Luer Lock	1,4 x 15 17G x 5/8" *
MPLL1420	Luer Lock	1,4 x 20 17G x 3/4" *
MPLL1425	Luer Lock	1,4 x 25 17G x 1" *
MPLL1435	Luer Lock	1,4 x 35 17G x 1" 3/8" *
MPLL1610	Luer Lock	1,6 x 10 16G x 3/8" *
MPLL1615	Luer Lock	1,6 x 15 16G x 5/8" *
MPLL1620	Luer Lock	1,6 x 20 16G x 3/4" *
MPLL1625	Luer Lock	1,6 x 25 16G x 1" *
MPLL1635	Luer Lock	1,6 x 35 16G x 1" 3/8" *
MPLL1815	Luer Lock	1,8 x 15 15G x 5/8" *
MPLL1820	Luer Lock	1,8 x 20 15G x 3/4" *
MPLL1825	Luer Lock	1,8 x 25 15G x 1" *
MPLL1835	Luer Lock	1,8 x 35 15G x 1" 3/8" *
MPLL1830	Luer Lock	1,8 x 30 15G x 1" 1/4" *
MPLL2020	Luer Lock	2,0 x 20 14G x 3/4" *
MPLL2025	Luer Lock	2,0 x 25 14G x 1" *
MPLL2030	Luer Lock	2,0 x 30 14G x 1" 1/4" *
MPLL2035	Luer Lock	2,0 x 35 14G x 1" 3/8" *
MPLL2040	Luer Lock	2,0 x 40 14G x 1" 5/8" *

* 100 pcs. package

MEDIOne



MEDI ONE NEEDLE (ALUMINIUM)

Code	Lock	Sizes
M1LL1213	Luer Lock	1,2 x 13 18G x 5/8" *
M1LL1220	Luer Lock	1,2 x 20 18G x 3/4" *
M1LL1225	Luer Lock	1,2 x 25 18G x 1" *
M1LL1240	Luer Lock	1,2 x 40 18G x 1" 1/2" *
M1LL1613	Luer Lock	1,6 x 13 16G x 5/8" *
M1LL1620	Luer Lock	1,6 x 20 16G x 3/4" *
M1LL1625	Luer Lock	1,6 x 25 16G x 1" *
M1LL1630	Luer Lock	1,6 x 30 16G x 1/4" *
M1LL1635	Luer Lock	1,6 x 35 16G x 1" 3/8" *
M1LL1640	Luer Lock	1,6 x 40 16G x 1" 1/2" *
M1LL2013	Luer Lock	2,0 x 13 14G x 5/8" *
M1LL2020	Luer Lock	2,0 x 20 14G x 3/4" *
M1LL2025	Luer Lock	2,0 x 25 14G x 1" *
M1LL2035	Luer Lock	2,0 x 35 14G x 1" 3/8" *
M1LL2040	Luer Lock	2,0 x 40 14G x 1" 1/2" *

MEDIOne Plastic



MEDI ONE NEEDLE (PLASTIC)

Code	Lock	Sizes
MM25058	Luer Lock	0,50 x 17 25 G x 5/8" *
MM23058	Luer Lock	0,60 x 17 23 G x 5/8" *
MM23100	Luer Lock	0,60 x 25 23 G x 1" *
MM23114	Luer Lock	0,60 x 32 23 G x 1" 1/4" *
MM22114	Luer Lock	0,70 x 32 22 G x 1" 1/4" *
MM21058	Luer Lock	0,80 x 17 21 G x 5/8" *
MM21100	Luer Lock	0,80 x 18 21 G x 1" *
MM21112	Luer Lock	0,80 x 19 21 G x 1" 1/2" *
MM21200	Luer Lock	0,80 x 20 21 G x 2" *
MM20100	Luer Lock	0,90 x 25 20 G x 1" *
MM19058	Luer Lock	1,0 x 16 19 G x 5/8" *
MM19100	Luer Lock	1,0 x 25 19 G x 1" *
MM19112	Luer Lock	1,0 x 38 19 G x 1" 1/2" *
MM19200	Luer Lock	1,0 x 51 19 G x 2" *
MM18112	Luer Lock	1,2 x 38 18 G x 1" 1/2" *
MM16100	Luer Lock	1,6 x 25 16 G x 1" *
MM16112	Luer Lock	1,6 x 38 16 G x 1" 1/2" *

* 100 pcs. package

