

## First class semen delivery bag



### NEW SILVER ROLL SACCO BAG

#### Silver SACCO bag

The aluminium version provides greater protection to the seminal fluid and ensures a top-quality image in terms of semen packaging.

#### Transparent SACCO bag

The transparent SACCO bag can be customized and is better in terms of quality of the plastic layers and cost-effectiveness.



Available also in transparent version

The SACCO bag revolutionised the AI industry and is the most effective genetic transfer method currently available.

#### Advantages of the new SACCO bag

- Designed to synchronise with the sows' natural physiological uterine contractions, thus avoiding any forced intake.
- Increases the semen dose shelf life thanks to the excellent surface-to-volume ratio, which reduces sperm cell degradation.
- Unique combination of two medical-grade, sperm-friendly, plastic layers (available both in aluminium and transparent material).
- Designed to be easily processed, labelled, transported and stored.
- More doses can be packed, loaded onto shipping containers and then stored in farms. Lower costs.
- Aseptically manufactured bag.
- Economically packaged in 500-piece rolls, resulting in lower shipping and storage costs.
- Possibility to customize the SACCO bag graphics, in order to show a better boar stud image to the final customer.



Re-closable tip with easier opening: a unique feature of the new SACCO bag



SACCO system

Re-closable tip allowing for semen analysis even before insemination.

User-friendly tip opening, easily connectable to the Medi Chimica catheters.

The SACCO bag is the most cutting-edge semen packaging system.

The SACCO range now includes two models, SACCO rolls and SACCO+, both available in aluminium as well as in transparent material.

#### SACCO

Code	Description
FM1000	SACCO bag silver roll
FM1001	SACCO bag transparent roll
FM1002	SACCO bag silver single
FM1003	SACCO bag transparent single